



**NEW**

# TableTop PZT™

## Compact Hard-Mount Vibration Cancellation System





*Sub-Hz compact active hard-mount vibration cancellation system for small precision instruments incorporates TMC's patented Stacis® technology*

### Features & Benefits

- ▶ **Active hard-mount, no soft air suspension**
- ▶ **Incorporates patented Stacis® technology**
- ▶ **Sub-Hz vibration cancellation, vertical and horizontal**
- ▶ **Lightweight, compact design**
- ▶ **Ideal for small precision instruments**
- ▶ **Extended stroke piezoelectric actuators, up to 60 µm**
- ▶ **Simple, robust, & cost-effective**
- ▶ **Payload capacity: up to 300 lb.**

**TableTop PZT™ Transmissibility**



### Ordering Chart

Catalog Number	Description	Pricing
25-406	TableTop PZT, 21.5 x 25.0 x 3.63 inches	Contact TMC

# TMC Products AT-A-GLANCE

## Laboratory Tables & TableTop™ Platforms



## Optical Tops, Breadboards, & Supports



## STACIS® Active Piezoelectric Vibration Cancellation System



## Floor Platforms for Raised & Solid Floors




## Electro-Damp® Active Vibration Cancellation System



## Pneumatic Vibration Isolators for OEM Applications



## Magnetic Field Cancellation & Electric Field Shielding Systems



## Acoustic Enclosures & Precision Structures



[www.techmfg.com](http://www.techmfg.com)

Copyright © 2007 by Technical Manufacturing Corporation. All rights reserved. All registered marks and trademarks are property of their respective holders.